



## ROLLZ FLEX - THE INCOGNITO WALKING FRAME



**SKU:** 1510-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Rollz Flex is a world first for its unique design and function and is perfect for those who need some support of a walking frame but want an alternative look to a traditional walker. Its large bag can safely carry up to 20kg so is a great companion for any mission. Light, strong, easy to steer and stylish.

### FEATURES

- Walking with large bag - the Flex boasts a large 25 litre shopping bag. The bag's centrally located position on the frame means that even when loaded with up to 20kg of shopping the frame will not tip. Compare this to other frames with the bag or basket located on the front of the frame which can be unsafe. This bag also includes an inner bag which can be easily lifted in and out.
- The highest quality fittings including sustainable, hard-wearing waterproof material on the charcoal bag.
- Comfortable padded and stylish seat using faux leather in black.
- Soft tyres made from a PU/foam mix for natural suspension
- 8inch wheels and short wheel base which make it very manoeuvrable and kind to your back

### ATTACHMENTS

[Click here](#) to view our full range of Walker/Frame attachments.

<https://youtu.be/iyW2G1ooZQ0>

## SPECIFICATIONS

- Rollz Flex: 7.7kg
- Rollz Flex bag: 25 litre
- Maximum carrying load: 125 kg
- Minimum and maximum height of user: 1.55- 1.95cm
- Height: 80-100cm
- Length: 60 cm
- Width: 62.3 cm
- Folded width: 33cm
- Seat height: 62 cm
- Seat width: 44 cm
- Diameter wheel: 8 inches

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)

