

ROLLZ MOTION ELECTRIC WALKER/WHEELCHAIR



SKU: 13101

PRODUCT DESCRIPTION

Purpose & Benefits of the Rollz Motion Electric

The Rollz Motion Electric is a revolutionary mobility aid that seamlessly combines a rollator and wheelchair, now with the added convenience of electric propulsion. Designed for individuals who need versatile support, this device allows users to walk comfortably with a rollator and switch effortlessly to a wheelchair when needed. With its electric motor assistance, the Rollz Motion Electric makes uphill and long-distance travel easier, offering enhanced mobility, independence, and freedom. Ideal for both indoor and outdoor use, it's perfect for those who want a compact, versatile solution for varied mobility needs.

Key Features

- **Dual Functionality:** Combines a rollator and wheelchair in one device, allowing users to switch between walking support and seated transport with ease.
- **Electric Propulsion Assistance:** Equipped with an electric motor to assist during long walks or uphill climbs, reducing the physical strain on the user.
- **Compact and Foldable Design:** The Rollz Motion Electric folds easily for transport and storage, making it perfect for travel or when space is limited.
- **Comfortable Ergonomic Seat:** Features a padded seat and adjustable backrest, providing comfort during extended use in the wheelchair mode.
- **Height Adjustable Handles:** The rollator's handles are adjustable to suit different user heights, ensuring proper posture and comfortable use.
- **All-Terrain Wheels:** Large wheels provide smooth navigation over various surfaces, making it suitable for indoor and outdoor use.
- **Battery-Powered:** Long-lasting battery life ensures extended use of the electric assistance, with easy charging and quick battery swaps when needed.

Specifications

- **Weight Capacity:** 125 kg
- **Weight (Without Battery):** 23.5 kg
- **Battery Life:** Up to 15 km of electric-assisted use on a single charge
- **Seat Width:** 45 cm
- **Wheel Size:** Large front and rear wheels for all-terrain capability
- **Folded Dimensions:** Compact and easy to store
- **Cleaning:** Wipe clean with a damp cloth
- **Usage:** Ideal for individuals who require both walking assistance and a wheelchair option, suitable for both home and outdoor use

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au