


ROLLZ MOTION RHYTHM WALKER



SKU: 16101

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Rollz Motion Rhythm is a walker with three cues for people with Parkinson's who experience freezing or a disrupted walking pattern. The laser line prompts the user to start a movement, while the metronome and the vibration in the handles support the walking rhythm. To start a cue, simply press the button on the handle.

FEATURES

- The cues on this device are highly customizable, allowing the user to select and adjust one, two, or all three to best meet their individual needs. Their rhythm and pattern can also be modified. You have the option to set the cues to run for a certain number of beats (8, 16, 32) or continuously.
- The audio and tactile cues can be discretely activated via your smartphone, with the audio potentially being streamed to earplugs through Bluetooth.
- The device also features a highly visible, bright green laser line.
- Its large wheels and relatively shorter wheelbase offer easy steering control.
- Ergonomically designed inverted handles lessen the strain on fingers and wrist joints.
- The device can be easily collapsed for convenient transportation.

SPECIFICATIONS

Rollator weight: 11.6 kg

Wheelchair weight: 3.4 kg

Maximum weight capacity: 125 kg
User height: 1.60 – 1.90 m
Seat height: 55 cm
Handle height range: 87 – 98 cm
Seat size rollator: 20 x 45 cm
Seat size wheelchair: 42 x 45 cm
Overall width: 67 cm
Overall length: 63 cm
Width folded: 31 cm
Front-wheel diameter: 8 inch
Rear-wheel diameter: 12 inch

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au