

## ROLLZ MOTION RHYTHM WALKER

**\$3,089.00**



**SKU: 16101**



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Rollz Motion Rhythm is a walker with three cues for people with Parkinson's who experience freezing or a disrupted walking pattern. The laser line prompts the user to start a movement, while the metronome and the vibration in the handles support the walking rhythm. To start a cue, simply press the button on the handle.

### FEATURES

- The cues on this device are highly customizable, allowing the user to select and adjust one, two, or all three to best meet their individual needs. Their rhythm and pattern can also be modified. You have the option to set the cues to run for a certain number of beats (8, 16, 32) or continuously.
- The audio and tactile cues can be discretely activated via your smartphone, with the audio potentially being streamed to earplugs through Bluetooth.
- The device also features a highly visible, bright green laser line.
- Its large wheels and relatively shorter wheelbase offer easy steering control.
- Ergonomically designed inverted handles lessen the strain on fingers and wrist joints.
- The device can be easily collapsed for convenient transportation.

### SPECIFICATIONS

Rollator weight: 11.6 kg

Wheelchair weight: 3.4 kg

Maximum weight capacity: 125 kg

User height: 1.60 – 1.90 m

Seat height: 55 cm

Handle height range: 87 – 98 cm

Seat size rollator: 20 x 45 cm

Seat size wheelchair: 42 x 45 cm

Overall width: 67 cm

Overall length: 63 cm

Width folded: 31 cm

Front-wheel diameter: 8 inch

Rear-wheel diameter: 12 inch