



ROLLZ MOTION - SLOW DOWN BRAKE

\$260.00



SKU: 10108

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The drag brakes offer continuous resistance to allow the user to control the speed of the rollator.

These slow down brakes are mounted above the wheels on both sides of the frame.

The desired resistance can be adjusted by tightening the adjustment screw and can be set to zero when used in wheelchair mode.

The drag brakes are compatible with any Rollz Motion rollator but it is not recommended to be used with pneumatic tyres.

The drag brakes come in a set of two.

Interested in the Rollz Performance 2 in 1 Walker/Wheelchair - [click here](#).