




ROLLZ MOTION - SLOW DOWN BRAKE



SKU: 10108

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The drag brakes offer continuous resistance to allow the user to control the speed of the rollator.

These slow down brakes are mounted above the wheels on both sides of the frame.

The desired resistance can be adjusted by tightening the adjustment screw and can be set to zero when used in wheelchair mode.

The drag brakes are compatible with any Rollz Motion rollator but it is not recommended to be used with pneumatic tyres.

The drag brakes come in a set of two.

Interested in the Rollz Performance 2 in 1 Walker/Wheelchair - [click here](#).

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au