



## ROLLZ RHYTHM PARKINSON ROLLATOR WALKER

**\$3,296.00**



**SKU: 16102**

## PRODUCT DESCRIPTION

### Purpose & Benefits

The Rollz Rhythm Parkinson Rollator Walker is specially designed for people with Parkinson's disease and related gait disorders. This advanced mobility aid features three integrated cueing modes—laser line projection, vibration in the handles, and an audible metronome—to help initiate movement and reduce freezing episodes during walking.

Its reverse braking system adds an extra layer of safety and control, while the ergonomic frame and padded seat provide comfort during rest. With adjustable handle height, supportive wheels, and a compact foldable frame, the Rollz Rhythm is an ideal rollator for use at home or in care settings where stability, safety, and neurological support are essential.

### Features

- Designed for users with Parkinson's or neurological gait disorders
- Built-in cueing: laser, metronome, and handle vibration
- Adjustable cue timing, frequency, and intensity
- Reverse braking system (squeeze to move) for added control

- Ergonomic, height-adjustable handles
- Padded seat and integrated bag
- Compact folding design for transport and storage
- Ideal for indoor and light outdoor use

## **Specifications**

- Maximum User Weight: 125 kg
- Overall Width: 67 cm
- Overall Length: 63 cm
- Width Folded: 31 cm
- Rollator Weight: 11.6 kg
- Seat Height: 55 cm
- Handle Height Range: 87-98 cm
- Suitable User Height: 1.60 m – 1.90 m
- Seat Size (Rollator Mode): 20 × 45 cm
- Seat Size (Wheelchair Mode): 42 × 45 cm
- Front Wheel Diameter: 8 inch
- Rear Wheel Diameter: 12 inch