

ROLLZ RHYTHM PARKINSON ROLLATOR WALKER



SKU: 16102

PRODUCT DESCRIPTION

Purpose & Benefits

The Rollz Rhythm Parkinson Rollator Walker is specially designed for people with Parkinson's disease and related gait disorders. This advanced mobility aid features three integrated cueing modes—laser line projection, vibration in the handles, and an audible metronome—to help initiate movement and reduce freezing episodes during walking.

Its reverse braking system adds an extra layer of safety and control, while the ergonomic frame and padded seat provide comfort during rest. With adjustable handle height, supportive wheels, and a compact foldable frame, the Rollz Rhythm is an ideal rollator for use at home or in care settings where stability, safety, and neurological support are essential.

Features

- Designed for users with Parkinson's or neurological gait disorders
- Built-in cueing: laser, metronome, and handle vibration
- Adjustable cue timing, frequency, and intensity
- Reverse braking system (squeeze to move) for added control

- Ergonomic, height-adjustable handles
- Padded seat and integrated bag
- Compact folding design for transport and storage
- Ideal for indoor and light outdoor use

Specifications

- Maximum User Weight: 125 kg
- Overall Width: 67 cm
- Overall Length: 63 cm
- Width Folded: 31 cm
- Rollator Weight: 11.6 kg
- Seat Height: 55 cm
- Handle Height Range: 87-98 cm
- Suitable User Height: 1.60 m – 1.90 m
- Seat Size (Rollator Mode): 20 × 45 cm
- Seat Size (Wheelchair Mode): 42 × 45 cm
- Front Wheel Diameter: 8 inch
- Rear Wheel Diameter: 12 inch

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au