



ROTHCARE 3 WHEEL ROLLATOR

\$166.00



SKU: RC-SC5041E

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Rothcare 3-Wheel Rollator is a lightweight and portable rollator that is ideal for people who need a little extra support while walking. It features a durable aluminium frame and weighs only 5kg, making it easy to lift and transport.

FEATURES

- Adjustable height handlebars
- A basket for carrying items
- Lockable brakes
- Large wheels for stability on uneven surfaces

The Rothcare 3-Wheel Rollator is ideal for people of all ages and abilities who need a little extra support while walking. It is especially well-suited for people who have difficulty with balance or who need a rollator to improve their mobility.

SPECIFICATIONS

- Large Wheels 200 x 50 mm solid
- Overall width 65 cm unfolded
- Folded width 30cm
- Length 59cm
- Handle height range 82 – 97 cm
- Weight of basket & Tray 700 grams
- Max user weight 110 kg

