




ROTHCARE LIGHTWEIGHT ALUMINIUM WALKER

\$145.00 - \$150.00 Price range: \$145.00 through \$150.00



SKU: RC-CM462-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Rothcare 4 Wheel Rollator available in both 6" (CM462) and 8" (CM463) tyre models, with height adjustable seat to 2 fixed positions. Easy to fold for storage and transport.

FEATURES

- Lightweight aluminium frame
- Durable 8" (CM463) or 6" (CM462) wheels
puncture free foam tyres
- Easy to fold for storage & transportation
- Padded seat with removable folding backrest
- Cane holder
- 2 fixed seat height positions

SPECIFICATIONS

	CM462 6" Red	CM463 8" Blue
Seat Width	30cm	30cm
Seat Depth	30cm	30cm
Seat Height	47-53cm	50-56cm
Width at Base	60cm	62cm
Depth at Base	68cm	70cm

	CM462 6" Red	CM463 8" Blue
Height	77-90cm	79-92cm
Weight Capacity	115kg	115kg
Total Weight	5.8kg	6.3kg