



ROTHCARE BOSTON SCOOTER



SKU: RC-SKPL00113-G-H

PRODUCT DESCRIPTION

Purpose & Benefits of Rothcare Boston Scooter

The **Rothcare Boston Scooter** is a stylish and compact mobility scooter designed to offer comfort, control, and convenience for everyday use. With intuitive controls and a practical design, the Boston is ideal for users seeking an easy-to-operate scooter for both indoor and outdoor environments.

Engineered for stability and ease of manoeuvrability, this scooter features a padded swivel seat, adjustable armrests, and a tiller that adapts to suit individual needs. The Boston's compact size makes it perfect for tight spaces, shopping centres, and public transport, while its smooth ride and reliable performance support users in maintaining independence.

Key Features:

- Compact 4-wheel design for stable and easy handling
- Swivel seat with padded support for comfort and access
- Adjustable tiller and armrests for a tailored fit
- Front basket for convenient storage on the go
- Suitable for indoor and outdoor use



Ideal for everyday errands, short trips, and tight spaces

Specifications:

Overall Length 1020 mm Overall Width 610 mm Overall Height 900 mm **Turning Radius** 1200 mm 10 km Max range KM Max Speed 8 km/h **Ground Clearance** 80 mm 24v/12 Ah **Battery** Motor Size 220w hp Charger 2Amp A 48 kg **Product Weight Total Weight Without Battery** 38 kg Max User Weight 136 kg **Brand** Rothcare

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au



NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au