



ROTHCARE BOSTON SCOOTER

\$1,951.03



SKU: RC-SKPL00113-G-H

PRODUCT DESCRIPTION

Purpose & Benefits of Rothcare Boston Scooter

The **Rothcare Boston Scooter** is a stylish and compact mobility scooter designed to offer comfort, control, and convenience for everyday use. With intuitive controls and a practical design, the Boston is ideal for users seeking an easy-to-operate scooter for both indoor and outdoor environments.

Engineered for stability and ease of manoeuvrability, this scooter features a padded swivel seat, adjustable armrests, and a tiller that adapts to suit individual needs. The Boston's compact size makes it perfect for tight spaces, shopping centres, and public transport, while its smooth ride and reliable performance support users in maintaining independence.

Key Features:

- Compact 4-wheel design for stable and easy handling
- Swivel seat with padded support for comfort and access
- Adjustable tiller and armrests for a tailored fit
- Front basket for convenient storage on the go
- Suitable for indoor and outdoor use

- Ideal for everyday errands, short trips, and tight spaces

Specifications:

Overall Length	1020 mm
Overall Width	610 mm
Overall Height	900 mm
Turning Radius	1200 mm
Max range KM	10 km
Max Speed	8 km/h
Ground Clearance	80 mm
Battery	24v/12 Ah
Motor Size	220w hp
Charger	2Amp A
Product Weight	48 kg
Total Weight Without Battery	38 kg
Max User Weight	136 kg
Brand	Rothcare