




ROTHCARE LIGHTWEIGHT CARBON FIBRE ROLLATOR

\$850.00



SKU: RC-INNU-W2301

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The RothCare super lightweight carbon fibre frame rollator has a cross folding frame for ease of loading into a car or storage with large wheels for outdoor or indoor use and stands when folded. It includes a removable handy storage bag, a removable padded backrest and a comfortable seat plus soft touch brakes.

FEATURES

- Folding frame for transportation with locking catch.
- Large wheels.
- Removable handy storage bag.
- Removable padded backrest.
- Comfortable seat.
- Soft touch brakes.

SPECIFICATIONS

- Width inside hand-grips 46cm
- Overall width 58cm
- Overall depth 70cm
- Seat Area 43 x 20 cm
- Seat Height 60cm

- Adjustable handle height 82 – 94 cm
- Front wheel 8 inch (PVC)
- Rear wheel 8 inch (PVC)
- Max. user weight 120kg
- Weight only 5.8 kgs including the bag