



ROTHCARE LITE RYDER CARBON FIBRE FOLDING POWERCHAIR



SKU: RC-INNU-N5909

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Rothcare Lite Rider Carbon Fibre Folding Power Chair, a remarkably lightweight mobility solution that weighs a mere 13kg, inclusive of the battery. This cutting-edge power wheelchair is designed to offer convenience and accessibility to individuals who may typically struggle with lifting heavier mobility devices into their car boot. In fact, it is even lighter than many non-motorised manual wheelchairs available in the market.

Crafted with a frame constructed from ultra-lightweight carbon fibre, this power wheelchair harnesses the renowned properties of this material, known for its exceptional lightness. The result is a power chair that opens up a world of possibilities, particularly in flat terrains such as shopping centres, hotels, airports, and galleries. However, it is important to note that this chair is not intended for navigating steep gradients or rough terrain.

Experience the freedom and ease of the Rothcare Lightweight Folding Power Wheelchair—a reliable and lightweight choice for effortless mobility in urban environments.

FEATURES

- Features a super lightweight carbon fibre frame for enhanced portability.
- Folds effortlessly in a single step, making it convenient for transportation.
- Equipped with a comfortable padded seat and backrest for added comfort during use.
- Offers a folding backrest for easy storage and compactness.
- Includes a lightweight removable Lithium battery for a convenient power supply.
- Features a compact and user-friendly joystick for easy maneuverability.
- The entire chair, including the battery, weighs less than 13kg, ensuring lightweight mobility.

- Comes with a USB fitting in the joystick and a removable battery pack for added convenience.
- The battery pack is fitted with a charging port, allowing for off-board charging.
- Provides an extra mounting point for an optional additional battery pack, offering extended usage time.

SPECIFICATIONS

- Overall width: 53cm
- Overall length: 90cm
- Overall height: 86cm
- Seat width: 43cm
- Seat depth: 40cm
- Seat height from the ground: 46cm
- Backrest height: 40cm
- Speed: 4.5kph
- Turning radius: 90cm
- Gradient 3 degrees: (1:19)
- Range: 10km
- Battery: 24V / 10Ah Lithium-ion
- Motor DC: 24V, 150W brushless DC
- Controller: Shanghai Micon
- Brake: Electromagnetic
- Front-wheel: 6"
- Rear-wheel: 10.6"
- Safe Working Load: 100kg
- Weight: w/out battery 11.85kg
- Complete Overall Weight with battery: 13kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au