



ROTHCARE LIGHTWEIGHT SELF PROPELLED WHEELCHAIR

\$638.00 - \$686.00Price range: \$638.00 through \$686.00



SKU: RC-AL-001JSP-16-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Rothcare Lightweight Self-Propelled Wheelchair is a stylish and lightweight wheelchair that is ideal for people who want to maintain their independence and mobility. It features a durable aluminium frame and weighs only 12.25kg, making it easy to lift and transport.

FEATURES

- · Lightweight and portable aluminium frame
- Weighs only 12.25kg
- Padded seat and backrest
- Full-length armrests with attendant brakes
- Swing-away footrests
- Removable leg rests
- Quick release 22" wheels
- Washable padded cushion set which is removable
- Modern ergonomically designed frame

SPECIFICATIONS

• Seat Width: 41cm (16.1"); 46cm (18.1"); 52cm (20.4")

• Length: 99cm

Overall width: 68cmOverall Height: 93cmSeat Depth: 43cm

• Seat Height from ground: 49cm



Backrest Height: 46cmWeight Capacity: 125kgTotal Weight: 12.25kg