



## **ROTHCARE LIGHTWEIGHT SELF PROPELLED WHEELCHAIR**

**\$708.95 - \$739.50** Price range: \$708.95 through \$739.50



**SKU:** RC-AL-001JSP-16-H

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Rothcare Lightweight Self-Propelled Wheelchair is a stylish and lightweight wheelchair that is ideal for people who want to maintain their independence and mobility. It features a durable aluminium frame and weighs only 12.25kg, making it easy to lift and transport.

### **FEATURES**

- Lightweight and portable aluminium frame
- Weighs only 12.25kg
- Padded seat and backrest
- Full-length armrests with attendant brakes
- Swing-away footrests
- Removable leg rests
- Quick release 22" wheels
- Washable padded cushion set which is removable
- Modern ergonomically designed frame

### **SPECIFICATIONS**

- Seat Width: 41cm (16.1"); 46cm (18.1"); 52cm (20.4")
- Length: 99cm
- Overall width: 68cm
- Overall Height: 93cm
- Seat Depth: 43cm
- Seat Height from ground: 49cm

- Backrest Height: 46cm
- Weight Capacity: 125kg
- Total Weight: 12.25kg