




## ROTHCARE LIGHTWEIGHT TRANSIT WHEELCHAIR

**\$574.00 - \$606.00** Price range: \$574.00 through \$606.00



**SKU:** RC-AL-001JTR-16-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Rothcare Lightweight Transit Wheelchair is a sturdy and lightweight aluminium wheelchair with a full range of individual adjustments. It is the perfect wheelchair for people who need a wheelchair that is both portable and customisable.

### FEATURES

- Sturdy and lightweight aluminium frame
- Full range of individual adjustments, including seat base depth and height, backrest sling height and contouring, push handle height, armrest depth and height, and footplate height and angle
- Attendant brakes
- Quick release wheels
- Anti-tip wheels
- Solid tyres front and rear
- Optional stump-support
- Available in 3 different sizes: 40.6cm (16") seat (Grey frame), 45.72cm (18") seat (Blue frame), and 50.8cm (20") seat (Red frame)

### SPECIFICATIONS

Parameter	RC-AL-001JTR-16	RC-AL-001JTR-18	RC-AL-001JTR-20
Seat width	40.6cm (16")	45.72cm (18")	50.8cm (20")
Overall width	55cm	60cm	65cm
Seat depth	40-44cm	40-44cm	40-44cm
Overall height	95-103cm	95-103cm	95-103cm

<b>Parameter</b>	<b>RC-AL-001JTR-16</b>	<b>RC-AL-001JTR-18</b>	<b>RC-AL-001JTR-20</b>
Seat to floor	51-55cm	51-55cm	51-55cm
Weight	14.7(12.9)kg	14.8(13)kg	14.9(13.1)kg
Weight capacity	130kg	130kg	150kg