



# ROTHCARE LIGHTWEIGHT TRANSIT WHEELCHAIR

\$517.00 - \$545.00Price range: \$517.00 through \$545.00



SKU: RC-AL-001JTR-16-H

### PRODUCT DESCRIPTION

# **PURPOSE & BENEFITS**

The Rothcare Lightweight Transit Wheelchair is a sturdy and lightweight aluminium wheelchair with a full range of individual adjustments. It is the perfect wheelchair for people who need a wheelchair that is both portable and customisable.

### **FEATURES**

- Sturdy and lightweight aluminium frame
- Full range of individual adjustments, including seat base depth and height, backrest sling height and contouring, push handle height, armrest depth and height, and footplate height and angle
- Attendant brakes
- Quick release wheels
- Anti-tip wheels
- Solid tyres front and rear
- Optional stump-support
- Available in 3 different sizes: 40.6cm (16") seat (Grey frame), 45.72cm (18") seat (Blue frame), and 50.8cm (20") seat (Red frame)

# **SPECIFICATIONS**

<b>Parameter</b>	RC-AL-001JTR-16	RC-AL-001JTR-18	RC-AL-001JTR-20
Seat width	40.6cm (16")	45.72cm (18")	50.8cm (20")
Overall width	55cm	60cm	65cm
Seat depth	40-44cm	40-44cm	40-44cm
Overall height	95-103cm	95-103cm	95-103cm



<b>Parameter</b>	RC-AL-001JTR-16	RC-AL-001JTR-18	RC-AL-001JTR-20
Seat to floor	51-55cm	51-55cm	51-55cm
Weight	14.7(12.9)kg	14.8(13)kg	14.9(13.1)kg
Weight capacity	130kg	130kg	150kg