



ROTHCARE MULTI-PURPOSE CANE



SKU: RC-BVPRSC0202A-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Rothcare Multi-Purpose Walking Cane is a versatile and reliable walking aid that can be used in a variety of situations. It features a wide base with four rubber feet for added stability, and a height-adjustable shaft to accommodate users of different heights.

The Multi-Purpose Walking Cane can be used as a standard walking cane, but it also has a number of other features that make it ideal for use in specific situations. For example, the wide base makes it ideal for use on uneven or slippery surfaces, while the height-adjustable shaft makes it easy to adjust the cane to the desired height for climbing stairs or standing up from a chair.

FEATURES

- Standard or Wide base with four rubber feet for added stability
- Height-adjustable shaft to accommodate users of different heights
- Soft foam handle for a comfortable grip

The Rothcare Multi-Purpose Walking Cane is ideal for people of all ages and abilities who need a little extra support while walking. It is especially well-suited for people who have difficulty balancing or who need to use a walker on a regular basis.

SPECIFICATIONS

- Lightweight and durable total weight of 750g
- Max user weight of 120kgs
- A Carry Bag is included

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au