



ROTHCARE QUAD CANE



SKU: RC-SC3001Q-BK

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Rothcare Quad Cane is a durable and stable walking aid that provides four points of contact with the ground for added stability. It is ideal for people who have difficulty balancing or who need extra support while walking.

The Quad Cane features a wide base with four rubber feet that help to prevent slipping and falling. It also has a height-adjustable shaft to accommodate users of different heights. The handle is made of soft foam for a comfortable grip.

The Quad Cane is easy to use. Simply place the cane in front of you and then step forward with your stronger leg. As you step, be sure to keep your back straight and your core engaged. The Quad Cane will help to support your weight and improve your balance.

FEATURES

- Wide base with four rubber feet for added stability
- Soft foam handle for a comfortable grip

SPECIFICATIONS

- Height from 74 -97 cm high
- Material: Powder coated, steel lower section, aluminium upper section
- 4 x 35mm diameter anti-slip ferrules (12mm shaft)
- Overall weight 800 grams
- Approx foot base 18 long x 13 cm wide
- Max user weight 100 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au