



## ROTHCARE STAIR ASSIST CANE



**SKU:** RC-SAC0202A

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Rothcare Stair Assist Cane is a specialized walking cane designed to help people safely ascend and descend stairs. It features a unique quad design with four rubber feet for added stability, and a height-adjustable shaft to accommodate users of different heights.

### FEATURES

- A wider base than a standard cane for increased stability
- A comfortable handgrip with two different heights to choose from
- Adjustable height with 12-section dual pin press buttons
- Use on hard surface steps or carpet

The Stair Assist Cane is easy to use. When going upstairs, simply place the cane on the next step and then step up with your stronger leg. When going downstairs, lead with the cane and then step down with your weaker leg.

### SPECIFICATIONS

- Max user weight 136 kg
- Total weight 750g

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)