



ROTHCARE STAIR ASSIST CANE

\$159.00



SKU: RC-SAC0202A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Rothcare Stair Assist Cane is a specialized walking cane designed to help people safely ascend and descend stairs. It features a unique quad design with four rubber feet for added stability, and a height-adjustable shaft to accommodate users of different heights.

FEATURES

- A wider base than a standard cane for increased stability
- A comfortable handgrip with two different heights to choose from
- Adjustable height with 12-section dual pin press buttons
- Use on hard surface steps or carpet

The Stair Assist Cane is easy to use. When going upstairs, simply place the cane on the next step and then step up with your stronger leg. When going downstairs, lead with the cane and then step down with your weaker leg.

SPECIFICATIONS

- Max user weight 136 kg
- Total weight 750g