




## ROTHCARE UNDERARM CRUTCH



**SKU:** RC-SC1002G-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Rothcare Underarm Crutch

The Rothcare Underarm Crutch is designed to provide optimal support and stability for individuals recovering from injuries or surgeries, or for those with long-term mobility challenges. These crutches are crafted to ensure a safe and comfortable mobility aid, reducing the risk of falls and facilitating recovery. The underarm design helps to distribute weight evenly, relieving pressure on the wrists and hands. Suitable for both short-term and long-term use, Rothcare Underarm Crutches are essential for anyone seeking a reliable and comfortable support system to aid in their mobility and independence.

### Features

- **Ergonomic Underarm Padding:** Offers comfortable support, reducing the risk of underarm discomfort and fatigue during extended use.
- **Adjustable Height:** Easily adjustable to fit the height of the user, ensuring optimal comfort and effectiveness in use.
- **Durable Construction:** Made with high-quality materials to withstand regular use while providing stable and reliable support.
- **Non-Slip Rubber Tips:** Equipped with non-slip rubber tips for added safety, ensuring stability on various surfaces and preventing slips.
- **Lightweight Design:** Despite its durability, the crutch is lightweight, making it easier for users to move around without extra effort.
- **Easy to Clean:** The materials used allow for easy cleaning and maintenance, ensuring hygiene and longevity of the crutches.

## Specifications

- **Material:** High-strength aluminium for the frame, with rubber tips and padded underarm supports
- **Safe Working Load:** 100kg
- **Adjustable Height Range:** Medium: 157cm - 178cm; Large: 178cm - 198cm.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)