



ROYAL COMFORT WIDE HIGH BACK LUMBAR ADJUSTABLE SUPPORT CHAIR

\$1,073.00 Original price was: \$1,073.00.<u>\$966.00</u>Current price is: \$966.00.



SKU: AGE700W

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Royal Comfort Wide High Back Chair is suitable for those needing full back support. It is adjustable height with user weight capacity up to 350kg.

FEATURES

- Royal Comfort Healthcare Chair
- · High Back Plus Head Cushion
- Rear Castors
- Depth Adjustable
- Custom Designed Heavy-Duty Metal Frame
- Frames are Finished in a Durable Powder-coat
- Seat & Back Rests are Solid Ply Timber
- Foams are Dunlop Enduro Fire Retardant
- Padded Armrests
- Excellent Comfort & Support
- Commercial Grade Upholstery
- High-Quality Stitching
- Hand Built in Australia
- Depth Adjustment
- 350kg Weight Capacity
- 10-Year Warranty



COLOUR & STYLE

Some colours in the image have been discontinued, select a colour from the available list of options.

DIMENSIONS

Wide

Size W610 x D450mm (Adjustable Depth)

Floor to Seat Adjustment H440 - 525mm

Back Rest Height 750mm
Arm Rest Height 200mm