




ROYALE FLUSH PLUMBED TOILET/BIDET LIFT SEAT

\$5,767.00



SKU: ROYWHI

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Royale Flush Plumbed toilet/commode lift seat provides a reliable and secure solution for transferring to the toilet, enabling individuals to comfortably and safely use their restroom at home. With its Electronic Toilet seat riser, this innovative device effortlessly replicates the natural standing movement with a simple touch of a button, allowing users to easily transition to a standing position on their own. Additionally, the Royale Flush model is equipped with a built-in Bidet, offering an added level of comfort and luxury during use.

*Please note - this seat must be **plumbed** in for the best user experience*

FEATURES

- Safe toileting at home gives you independence
- Lowers you under power from a standing position until you are positioned over the toilet
- Raises you back up under power from a sitting position
- Safely use the toilet with reduced assistance or on your own
- Helps reduce manual handling risks for carers who provide toileting assistance
- Reduce the risk of falls and associated injuries, and reduce the likelihood of be unable to get up from the toilet
- Electric washing Bidet feature included (plumbing required)
- Remote control for ease of use
- Maximum Load weight 200kg

SPECIFICATIONS

Dimensions (L x W x H)

60.7 x 52.5 x 76.6cm

Weight	24kg
Max. Loading Weight	200kg
Distance between handrails	50.4cm
Support times for battery full charge	N/A Unit is plugged into power outlet
Working Voltage	110-240V adaptive voltage
Shape	33 degree up and down

SPECIFICATIONS

[Royale-Flush - Toilet Commode Seat Information](#)