




ROYALE MEDICAL HEAVY DUTY BATH TRANSFER BENCH



SKU: S11190

PRODUCT DESCRIPTION

Purpose & Benefits of Royale Medical Heavy Duty Bath Transfer Bench

The Royale Medical Heavy Duty Bath Transfer Bench is engineered to provide a safe, comfortable, and stable platform for individuals with mobility issues to transition into and out of the bathtub. This bench is especially beneficial for bariatric patients and those who find it difficult to stand in the shower, enhancing their ability to maintain hygiene independently in a secure manner.

Features

- **Sturdy and Rust-Resistant:** Made from aluminium, which is not only lightweight but also rust-resistant, ensuring durability and longevity.
- **Adjustable Height:** The legs of the bench can be adjusted to various heights, accommodating different user needs and bathtub heights.
- **High Weight Capacity:** Designed to safely support a substantial weight, suitable for bariatric patients.
- **Sliding and Swivel Seat:** The seat slides and rotates, making transfers smoother and less physically demanding for both users and caregivers.
- **Comfortable Seating:** Equipped with a padded seat, backrest, and armrests to ensure comfort during use.
- **Safety Features:** Includes a safety belt and large suction cups on the feet for additional stability and security.

Specifications

- **Safe Working Load:** 280 kg
- **Dimensions:**

- Overall Height: 82.5 cm - 97.5 cm
- Overall Length: 52 cm
- Overall Width: 78.5 cm
- Seat Height: 49 cm - 64 cm
- Seat Depth: 41 cm
- Seat Width: 68 cm
- Armrest Height: 27.5 cm
- Unit Weight: 6 kg

The Royale Medical Heavy Duty Bath Transfer Bench is an essential aid for individuals who require support while bathing. It provides a stable and secure environment, reducing the risk of slips and falls in the bathroom, thereby promoting independence and preserving dignity.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au