



SAEBOSTRETCH DYNAMIC HAND SPLINT



SKU: SBO-S1SS10SMR1-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The SaeboStretch Dynamic Hand Splint is a revolutionary device designed to help individuals with neurological impairments maintain or improve the functionality of their hands and wrists. Particularly beneficial for stroke survivors or those with conditions that cause spasticity, the SaeboStretch has a unique design that adapts to the tone and changing needs of the user, offering support and promoting joint alignment.

FEATURES

- **Dynamic Tension:** The SaeboStretch includes adjustable tension bands that gently guide the fingers into extension in response to involuntary flexion.
- **Customisable Fit:** The splint provides a customisable strapping system and a malleable frame that can be adjusted to accommodate various hand sizes and abnormalities.
- **Ergonomic Design:** With a focus on comfort and effectiveness, the splint is ergonomically designed to fit the contours of the hand and wrist, minimizing pressure points.
- **Soft and Breathable Cover:** The cover of the splint is soft to the touch, promoting comfort during prolonged use, and is made from a breathable material to help keep the skin dry.
- **Removable Components:** For easy cleaning and maintenance, the SaeboStretch features removable components, which can be hand-washed or replaced as needed.

The SaeboStretch Dynamic Hand Splint stands as a testament to Saebo's commitment to improving the lives of people with impaired hand function. This device is not only designed to assist with daily living but also to support rehabilitation and enable greater independence. By incorporating the latest in therapeutic design and materials, the SaeboStretch is a vital tool in the journey towards recovery and improved hand use.

SIZING INSTRUCTIONS

Measure across the proximal phalanges of digits 2-5 just proximal to the PIP joints. Please make sure

the fingers are in the adducted position. Do not wrap the tape down the lateral or medial side of the hand. Consider the length of the fingers when deciding between 2 different sizes. (i.e., shorter fingers = smaller size).

Sizing Chart

S	5.7cm - 7.0cm
M	7.0cm - 8.2cm
L	8.2cm - 10.0cm

Please ensure you measure correctly and get the right size as this item cannot be returned.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au