



SAEBOGLOVE COMPLETE KIT



SKU: SaeboGlove-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **SaeboGlove Kit** is a cutting edge, low-profile functional glove that assists with thumb and finger extension for clients with neurologic and orthopaedic impairments. It has an innovative design that features a spiral forearm splint that secures the wrist in a functional position. Various sized extension assist tensioners are located at the IP joints for all digits and assist with extending the fingers and thumb following grasping. The palm is exposed to increase breathability and ease for donning.

Features

- Low-profile functional glove that assists thumb and finger extension
- Spiral forearm splint secures wrist in a functional position
- Various-sized extension tensioners at the IP joints for all digits, allowing custom assist based on finger length and weakness
- Individual tensioners can be removed or adjusted to match user needs
- Palm is exposed to improve breathability and ease of donning
- Non-slip liner to reduce glove migration during use



Silicone-covered fingertips to improve grip during grasping activities

Specifications

- Sizes Available: Multiple sizes for both shell and liner; Left and Right-hand options
- Tensioners: Various sizes included, located at IP joints for fingers and thumb
- Wrist Splint: Spiral design to maintain functional wrist position
- Liner Material: Lycra material for expandability, non-slip liner
- Palm Exposure: Palm remains open (not covered) for breathability and ease of donning

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT



Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au