



SAEBOGLOVE COMPLETE KIT

\$1,103.00



SKU: SaeboGlove-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **SaeboGlove Kit** is a cutting edge, low-profile functional glove that assists with thumb and finger extension for clients with neurologic and orthopaedic impairments. It has an innovative design that features a spiral forearm splint that secures the wrist in a functional position. Various sized extension assist tensioners are located at the IP joints for all digits and assist with extending the fingers and thumb following grasping. The palm is exposed to increase breathability and ease for donning.

Features

- Low-profile functional glove that assists thumb and finger extension
- Spiral forearm splint secures wrist in a functional position
- Various-sized extension tensioners at the IP joints for all digits, allowing custom assist based on finger length and weakness
- Individual tensioners can be removed or adjusted to match user needs
- Palm is exposed to improve breathability and ease of donning
- Non-slip liner to reduce glove migration during use

- Silicone-covered fingertips to improve grip during grasping activities

Specifications

- **Sizes Available:** Multiple sizes for both shell and liner; Left and Right-hand options
- **Tensioners:** Various sizes included, located at IP joints for fingers and thumb
- **Wrist Splint:** Spiral design to maintain functional wrist position
- **Liner Material:** Lycra material for expandability, non-slip liner
- **Palm Exposure:** Palm remains open (not covered) for breathability and ease of donning