



SAEBOSTEP COMPLETE KIT V2

\$454.00



SKU: SBO-S3-KT-12-SD-N-1

PRODUCT DESCRIPTION

Purpose & Benefits

The **SaeboStep Complete Kit V2** is a lightweight drop foot brace designed to provide optimum foot clearance and support during walking. It replaces bulky, uncomfortable splints and less effective external braces, offering a more comfortable, durable, and adjustable solution for everyday mobility.

Features

- Strong, durable spectra cord quickly attaches to included eyelet hooks for effective foot lift
- BOA dial technology allows fast and precise adjustment of lift angle
- Hook-and-loop Velcro strapping system secures the brace to the ankle
- Convenient release function on dial for easy tension release throughout the day
- Compatible with most shoe styles, with optional accessory kits available for shoes without eyelets

Specifications

- **Included in Kit:** SaeboStep brace, 2 x eyelet kits (each includes 2 stainless steel C hooks, 2 x 3/16" screws, 2 x 1/4" screws, 2 x 5/16" screws, 2 x 11mm washer discs, 1 x Allen wrench)

- **Colour:** Black
- **Compatibility:** SaeboStep Eyelet Hook Kit and accessory kits for non-eyelet shoes