



SALJOL ALUMINIUM ROLLATOR



SKU: 42102-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Saljol Aluminium Rollator** is a versatile, lightweight mobility aid designed for comfort, safety and ease of use. Weighing only 6.4kg, it supports users up to 150kg and comes fully equipped with a backrest, storage bag, cane holder, bell and care cloth. With its strong aluminium frame, foldable design and premium reflectors, it's ideal for both indoor and outdoor use — delivering stability and confidence with every step.

FEATURES

- All-round performer crafted from durable, lightweight aluminium
- Includes backrest, storage bag, cane holder, bell and care cloth
- Optional one-handle braking system ideal for users with limited dexterity
- Folds easily for transport and storage
- Large front all-terrain tyres for smooth outdoor performance
- Wide handle height range up to 103cm for optimal ergonomic fit
- Reflective detailing enhances visibility and safety in low-light conditions
- Seat and handles are height adjustable for personalised comfort
- No assembly required ready to use straight out of the box

SPECIFICATIONS

| | | | | | Handle | | | |
|-------|--------|----------|-----------|-------|--------|---------|--------|----------|
| Model | Weight | User | Seat Seat | Seat | Height | Overall | Folded | Bag |
| | | Capacity | Heiaht | Width | | Width | Width | Capacity |
| | | 1 7 | - 5 | _ | Range | _ | _ | 1 , |



| Model | Weight | User Capacity | Seat Height | Seat Width | Handle Height Range | Overall Width | Folded Width | Bag Capacity |
|-------------------------------------|--------|------------------|----------------|---------------|---------------------------|------------------|-----------------|-----------------|
| Regular (AR62) | 6.4kg | 150kg | 62cm | 46cm | 81 – 103cm | 61cm | 21cm | 5kg |
| Small (AR54) | 6.3kg | 150kg | 54cm | 46cm | 76 - 88cm | 61cm | 21cm | 5kg |
| Wide (AR62W – coming 2025) | 6.5kg | 150kg | 62cm | 52cm | 81 – 103cm | 61cm | 21cm | 5kg |

MATERIALS

• Frame: Aluminium

• Handgrip: TPE with 20-30% cork

• Brake Handle: PA6 GF30 + TPE with 20-30% cork

Brake Pads: PA6 GF15
Seat: 100% polyester
Wheel Rim: PP GF15
Wheel Tyres: PU

Note: Accessories included. Specifications vary by model.

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom 343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777 elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771 reynella@medimart.com.au

VIC



Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au