

SALJOL ALUMINIUM ROLLATOR



SKU: 42102-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Saljol Aluminium Rollator** is a versatile, lightweight mobility aid designed for comfort, safety and ease of use. Weighing only 6.4kg, it supports users up to 150kg and comes fully equipped with a backrest, storage bag, cane holder, bell and care cloth. With its strong aluminium frame, foldable design and premium reflectors, it's ideal for both indoor and outdoor use — delivering stability and confidence with every step.

FEATURES

- All-round performer crafted from durable, lightweight aluminium
- Includes backrest, storage bag, cane holder, bell and care cloth
- Optional one-handle braking system — ideal for users with limited dexterity
- Folds easily for transport and storage
- Large front all-terrain tyres for smooth outdoor performance
- Wide handle height range up to 103cm for optimal ergonomic fit
- Reflective detailing enhances visibility and safety in low-light conditions
- Seat and handles are height adjustable for personalised comfort
- No assembly required — ready to use straight out of the box

SPECIFICATIONS

Model	Weight	User Capacity	Seat Height	Seat Width	Handle Height Range	Overall Width	Folded Width	Bag Capacity
-------	--------	---------------	-------------	------------	---------------------	---------------	--------------	--------------

Model	Weight	User Capacity	Seat Height	Seat Width	Handle Height Range	Overall Width	Folded Width	Bag Capacity
Regular (AR62)	6.4kg	150kg	62cm	46cm	81 – 103cm	61cm	21cm	5kg
Small (AR54)	6.3kg	150kg	54cm	46cm	76 – 88cm	61cm	21cm	5kg
Wide (AR62W – coming 2025)	6.5kg	150kg	62cm	52cm	81 – 103cm	61cm	21cm	5kg

MATERIALS

- **Frame:** Aluminium
- **Handgrip:** TPE with 20–30% cork
- **Brake Handle:** PA6 GF30 + TPE with 20–30% cork
- **Brake Pads:** PA6 GF15
- **Seat:** 100% polyester
- **Wheel Rim:** PP GF15
- **Wheel Tyres:** PU

Note: Accessories included. Specifications vary by model.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au