



## SALLYBOARD CURVED PATIENT TRANSFER BOARD

**\$105.00 - \$244.00** Price range: \$105.00 through \$244.00



**SKU:** MT SBC-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The SallyBoard Curved Patient Transfer Board is lightweight and allows easy patient transfers. The SallyBoard has a smooth top surface and curved edges.

A curved transfer board (or curved slide board) is a device used to help individuals with a disability or limited mobility transfer safely from one surface to another. The transfer board acts as a bridge between two surfaces, minimising the risk of slips and falls.

The SallyBoard Curved was specially developed to provide a light, chip and splinter-resistant alternative to other curved transfer boards.

Small and light, with a slight curve, it is popular not only in healthcare facilities but for home users – facilitating easy transfers from chair or wheelchair to car, toilet, bed, and more.

### FEATURES

- Constructed from lightweight HDE.
- Two squares underneath for grip and stability
- Good level of flexibility in the board where bed/barouche heights are not always at the same level
- Radiolucent: can be used for X-rays
- TGA registered
- Tested to Australian Standards AS/NZS ISO 10535:2011
- Straight option also available (SMMTSB)

## **SPECIFICATIONS**

- Height: 22.2cm
- Width: 77cm
- Thickness: 1cm
- Weight: 1.4kg
- SWL: 130kg when the gap is 20-30cm; 180kg when the gap is less than 20cm