

SAMMONS PRESTON FOAM WHEELCHAIR CUSHION



SKU: PAT-081

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Lightweight Cushion with Ties for Secure Fit - The Sammons Preston Wheelchair Foam Cushion is a lightweight cushion that has ties on the back to allow for a secure fit into the wheelchair. This cushion comes in multiple sizes to fit any wheelchair, providing individuals with a comfortable sitting surface while reducing discomfort.

High-Density Polyurethane Foam and Extra Strong Poly/Cotton Cover - The Sammons Preston Wheelchair Foam Cushion is made with high-density polyurethane foam and an extra strong poly/cotton cover that is both removable and washable. The high-density foam conforms to the shape of the body, providing a comfortable sitting surface while reducing pressure and discomfort. Measuring 16" deep, this cushion is an excellent investment for individuals who require additional support and comfort when sitting in their wheelchair.

Meets CA117 Standards - The Sammons Preston Wheelchair Foam Cushion meets CA117 standards, ensuring that it is a safe and reliable option for individuals who use wheelchairs. Its high-density foam and extra strong poly/cotton cover work together to reduce pressure and discomfort while sitting in a wheelchair. Plus, it is available in multiple sizes so it can fit any wheelchair.

FEATURES

- Lightweight cushion with ties for secure fit.
- High-density polyurethane foam and extra strong poly/cotton cover.
- Removable and washable cover meets CA117 standards.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au