




SAMMONS PRESTON HIP KIT III WITH 60CM HANDI REACHER

\$86.00



SKU: 081286467

PRODUCT DESCRIPTION

Purpose & Benefits

The **Sammons Preston Hip Kit III with 60 cm Handi Reacher** provides a complete set of aids to support recovery after hip replacement surgery. It is also ideal for seniors or individuals with arthritis who have difficulty bending. The kit promotes safe, independent dressing and bathing at home or in rehabilitation settings.

Features

- Sock Notch allows socks to be put on without bending
- Long shoe horn makes it easier to slide on shoes while protecting the backs
- Dressing stick helps with pulling up pants and adjusting clothing
- Long handled sponge allows easier washing of legs, feet and back
- 60 cm Handi Reacher reduces the need to bend past 90 degrees by retrieving items from the floor or hard-to-reach areas

Specifications

- **Reacher Length:** 60 cm

- **Kit Contents:** Handi Reacher, Sock Notch, Long Shoehorn, Dressing Stick, Long Handled Sponge