

SAVANAH TOILET SEAT RAISER

\$63.00 - \$90.00 Price range: \$63.00 through \$90.00



SKU: PAT-AA211-H

PRODUCT DESCRIPTION

SAVANAH TOILET SEAT RAISER PURPOSE & BENEFITS

Lightweight, all plastic construction, for strength and durability, the Savannah Toilet Seat Raiser has a gently contoured surface for extra comfort. The seat slopes downwards slightly from back to front, so as to allow easier lowering and raising, and has a deep splash guard to minimise soiling. Resistant to stains and odours, it is fully sealed and easy to clean.

- The large rear cut-out on the Savannah gives clearance for the coccyx and, combined with the gentle downward slope of the seat, makes rear personal hygiene possible.
- The wide front cut-out also allows easier access for personal cleansing.
- Available with a moulded plastic lid to create a more pleasing appearance.

SPECIFICATIONS

- Available in three height sizes, without or with Lid.
- Overall Length: 407mm (16")
- Overall Width: 370mm (14½") (excluding brackets)
- Aperture Length: 266mm (10½")
- Aperture Width: 207mm (8")
- Maximum User Weight: 190kg / 30st

You can also have a look at this toilet seat, [check here](#).