



## **SHEAR COMFORT DAY CHAIR WOOL OVERLAY WITH STRAPS**

**\$264.00 - \$299.00** Price range: \$264.00 through \$299.00



**SKU:** XDDCBS-H

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The **Shear Comfort XD1900 Day Chair Overlay with Straps** provides skin protection and superior comfort for individuals seated in day chairs or recliners for extended periods. Designed to reduce pressure, shear, and friction, this overlay helps maintain skin integrity while delivering a soft and supportive surface. The addition of securing straps ensures the overlay stays in place, offering stability and peace of mind during use.

### **Features**

- Constructed from Shear Comfort XD1900 high-density medical wool
- Helps redistribute pressure to prevent skin breakdown
- Reduces shear and friction for enhanced comfort
- Breathable, moisture-wicking, and temperature regulating
- Includes straps for secure placement on day chairs or recliners
- Machine washable for easy cleaning

## Specifications

- **Brand:** Healthcare Innovations Australia (Shear Comfort)
- **Product Type:** Day chair overlay with straps
- **Size:** 50cm x 140cm
- **Material:** XD1900 medical-grade wool
- **Colour:** Available in White or Blue