

SHEAR COMFORT HEEL PROTECTORS (PAIR)

\$116.00Online Only



SKU: XDHP-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Providing security and comfort to the heel and ankle areas, Shear Comfort Heel Protectors provide effective prevention of skin ulcers in at risk areas.

The wool wraps around the heel and ankle, while leaving the rest of the foot unencumbered, and is held in place by a Velcro strap. Seams are positioned to be away from any skin areas at risk of ulcers.

Sold in pairs, heel protectors are made out of XD1900 extra density material made of natural wool. XD1900 is manufactured in such a way as to improve the wool fibre density and therefore the protective properties of the heel protector.

Comes in a pair.