

## SHEAR COMFORT PRESSURE CARE BOOT





SKU: SBPBI -

## **PRODUCT DESCRIPTION**

# **Purpose & Benefits**

The **Shear Comfort Pressure Care Boot** is designed to provide protection, comfort, and pressure relief for people with delicate or at-risk skin. The boot's high cut and supportive structure help stabilise the foot and ankle, reducing shear and friction while keeping the skin safe. It is ideal for individuals with limited mobility, those prone to skin breakdown, or anyone requiring extra protection for foot health.

#### **Features**

- High-cut boot design for enhanced protection and support
- Lined with Shear Comfort XD1900 medical-grade wool for softness and durability
- Reduces shear and friction forces to safeguard fragile skin
- · Adjustable fastening straps allow for a secure and customised fit
- Breathable and moisture-wicking material supports skin integrity
- Machine washable for easy cleaning and care



# **Specifications**

• Brand: Healthcare Innovations Australia (Shear Comfort)

• Material: XD1900 medical-grade wool lining

• **Design:** Pressure care boot with high-cut style and adjustable straps

• Use: Protective footwear to reduce pressure, shear, and friction on vulnerable skin