



SHEAR COMFORT XD1900 CALF PLATE PROTECTORS (PAIR)

\$60.72



SKU: 2138

PRODUCT DESCRIPTION

Purpose & Benefits

The **Shear Comfort Calf Plate Protector** is specifically designed to cushion the often unforgiving calf plate regions of a wheelchair. Crafted from XD1900 high-density medical-grade wool, it offers superior pressure redistribution and maintains a healthy microclimate by efficiently managing heat and moisture. Elastic straps ensure the protectors stay securely in place, while easy machine washability and durability support long-term hygiene and use.

Features

- Made from XD1900 high-density wool (1,900 g/m²) for effective pressure redistribution
- Elastic straps for secure attachment to wheelchair calf plates
- Maintains a healthy skin microclimate—wool absorbs up to 30 % of its weight in moisture without feeling damp
- Machine washable up to 80 °C for hygienic care
- Guaranteed to withstand up to 50 washes with proper care

Specifications

- **Brand:** Shear Comfort (Healthcare Innovations Australia)
- **Product Type:** Calf Plate Protector (pair)
- **Material:** XD1900 medical-grade wool (fabric-backed)
- **Colour:** Blue
- **Washability:** Machine wash up to 80 °C, up to 50 washes guaranteed