



SHEAR COMFORT XD1900 CLASSIC OVERLAY - BLUE

\$275.00



SKU: XDCOF

PRODUCT DESCRIPTION

Purpose & Benefits

The **Shear Comfort Medical Sheepskin Overlay - Classic Size 75 cm** × **125 cm** provides plush, natural comfort and advanced pressure care when used on beds, chairs, or armchairs. Crafted from high-density XD1900 medical-grade wool, it offers superior pressure redistribution, a breathable and moisture-wicking microclimate, and the calming feel of pure wool for better rest and skin protection.

Features

- Made with XD1900 high-density medical-grade wool (1,900 g/m²), conforming to Australian Standard AS4480.1
- Soft, thick wool for enhanced comfort and pressure redistribution
- Maintains a healthy microclimate—absorbs up to 30 % of its weight in moisture vapour without feeling damp and allows heat and moisture to dissipate
- Promotes quality sleep—naturally advantageous over plastic-based solutions
- Machine washable up to 80 °C thanks to fabric-backed construction
- Guaranteed to withstand up to 50 washes when cared for correctly



Specifications

• Brand: Shear Comfort (Healthcare Innovations Australia)

• Product Type: Medical Sheepskin Overlay (Classic size)

• Material: XD1900 medical-grade wool (fabric-backed)

• Colour: Blue

• **Dimensions:** 75 cm × 125 cm

• Washability: Machine wash up to 80 °C; up to 50 washes guaranteed