




SHEAR COMFORT XD1900 PALM PROTECTOR

\$57.00Online Only



SKU: XDPALM

PRODUCT DESCRIPTION

Purpose & Benefits

The Shear Comfort XD1900 Palm Protector is designed to provide relief and protection for the palm and hand, especially for individuals who experience contractures, pressure points, or skin breakdown due to limited mobility or constant hand closure.

Crafted from high-density XD1900 medical-grade sheepskin, the palm protector cushions the hand and keeps the fingers comfortably open, helping prevent fingernail pressure against the palm and reducing the risk of pressure injuries.

Features

- Premium XD1900 sheepskin for high-density cushioning
- Helps prevent pressure sores and skin breakdown
- Supports hand contractures and protects sensitive palms
- Soft, moisture-wicking and breathable material
- Secured with adjustable Velcro strap
- Machine washable up to 80°C with Shear Comfort shampoo

- Hypoallergenic and gentle on fragile skin

Specifications

- Size: One size fits most
- Colour: Blue
- Material: XD1900 medical sheepskin (1900gsm)
- Washable: Up to 80°C