

## **SHEAR PRESSURE CARE ASSISTANT XD1900**



**SKU:** XDPCA



## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The **Shear Comfort Pressure Care Assistant** is designed to provide targeted relief and comfort for specific body areas—such as heels, hips, knees, shoulders, or head—while lying in bed. It securely attaches via linen flaps tucked under the mattress, keeping the wool protection exactly where it's needed. Made from high-density XD1900 medical-grade wool, it offers superior pressure redistribution, a healthy skin microclimate, and long-lasting durability.

### **Features**

- Linen flaps enable easy secure attachment under the mattress for precise placement
- Constructed with high-density XD1900 wool (1,900 g/m<sup>2</sup>) for effective pressure redistribution
- Promotes a healthy microclimate—wool absorbs up to 30 % of its weight in moisture vapor without feeling damp, allowing heat and moisture to dissipate
- Machine washable up to 80 °C thanks to fabric-backed construction for hygienic cleaning
- Guaranteed durability for up to 50 washes when cared for according to manufacturer's instructions

## Specifications

- **Brand:** Shear Comfort (Healthcare Innovations Australia)
- **Product Type:** Pressure Care Assistant (bed overlay)
- **Material:** XD1900 medical-grade wool with fabric backing
- **Colour:** Blue
- **Dimensions:** 70 cm x 30 cm
- **Washability:** Machine washable up to 80 °C; up to 50 washes guaranteed

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)