

DAYS SHOPPER 12 EXTRA WIDE TRANSIT WHEELCHAIR



SKU: MFI-VNW01AYA



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Shopper 12 Extra Wide 22" Wheelchair has all the same great features as the Standard Shopper Wheelchair, but also features a larger seat and removable armrests.

FEATURES

- Light 13.5kg frame (excluding leg-rests);
- Height adjustable;
- Larger seat;
- Removable armrests;
- Folding backrest;
- Attendant hand brakes;
- Rear wheel locks;
- Safety lap belt;
- Swing-away flip-up footplates; also
- Padded armrests.

SPECIFICATIONS

Weight of Chair:	13.5kg frame (15kg with legrests)
Overall Seat Width:	640mm
Seat width (Between Arms):	560mm
Seat Depth:	407mm

Backrest Height:	457mm
Legrest Height Range:	356-483mm (add cushion for additional height)
Overall Width (Footprint):	680mm
Overall Depth (Footprint):	940mm
Armrest Height (Seat to Top of Handle):	254mm
Push Handle Height:	960mm
Front Wheels:	203 x 25.4mm - Solid Castors
Rear Wheels:	305 x 57mm
Safe Working Load:	125kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au