



## **DAYS SHOPPER 12 EXTRA WIDE TRANSIT WHEELCHAIR**

**\$497.00**



**SKU: MFI-VNW01AYA**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Shopper 12 Extra Wide 22" Wheelchair has all the same great features as the Standard Shopper Wheelchair, but also features a larger seat and removable armrests.

### **FEATURES**

- Light 13.5kg frame (excluding leg-rests);
- Height adjustable;
- Larger seat;
- Removable armrests;
- Folding backrest;
- Attendant hand brakes;
- Rear wheel locks;
- Safety lap belt;
- Swing-away flip-up footplates; also
- Padded armrests.

### **SPECIFICATIONS**

Weight of Chair:	13.5kg frame (15kg with legrests)
Overall Seat Width:	640mm
Seat width (Between Arms):	560mm
Seat Depth:	407mm

Backrest Height:	457mm
Legrest Height Range:	356-483mm (add cushion for additional height)
Overall Width (Footprint):	680mm
Overall Depth (Footprint):	940mm
Armrest Height (Seat to Top of Handle):	254mm
Push Handle Height:	960mm
Front Wheels:	203 x 25.4mm - Solid Castors
Rear Wheels:	305 x 57mm
Safe Working Load:	125kg