



## SHOPRIDER ALLROUNDER MOBILITY SCOOTER

**\$4,818.00 - \$4,895.00** Price range: \$4,818.00 through \$4,895.00



**SKU:** 9LSP-H

## PRODUCT DESCRIPTION

### Features and Benefits

The Shoprider All Rounder is a high-performance mobility scooter built for comfort, control, and all-terrain reliability. With large 330mm wheels and a robust rear-wheel-drive motor, it provides smooth navigation across uneven surfaces while offering enhanced ground clearance.

Designed for user convenience, it features a looped steering tiller with intuitive controls, an electromagnetic brake system for safety, and a 1.5Hp motor that pairs with 40Ah batteries to deliver up to 40km range on a single charge. This scooter is ideal for users seeking a reliable, feature-rich solution for both everyday use and longer outings.

### Features

- Large 330mm wheels for superior ground clearance
- Rear-wheel drive and electromagnetic brake for control and safety
- 1.5Hp motor and 8Amp charger for optimal performance
- Comfortable seating and ergonomic looped tiller
- Travel range of up to 40km per charge

## Specifications

- **Length:** 1300mm
- **Width:** 640mm
- **Height:** 1200mm
- **Safe Climbing Angle:** 9°
- **Max Climbing Angle:** 12°
- **Maximum Load:** 159kg
- **Charger:** 8Amp
- **Motor Size:** 1.5Hp
- **Battery Capacity:** Minimum 40Ah
- **Maximum Range:** Up to 40km
- **Wheel Size:** 330mm
- **Turning Radius:** 1440mm
- **Weight (with battery):** 99kg