



## SHOPRIDER EXPLORER

**\$3,950.00**



**SKU: 888SE**

## PRODUCT DESCRIPTION

### Features and Benefits

The Shoprider Explorer is a refreshed and enhanced version of the trusted 888SE model, offering excellent comfort, control, and convenience for users seeking reliable personal mobility. The new model features a gas-assisted adjustable tiller, a high back contoured seat, and improved usability elements like a USB port and LED headlight.

It's designed with user experience in mind, incorporating adjustable armrests, a padded loop handlebar with fingertip controls, and multiple customisation points for height, angle, and positioning. The Explorer also includes an updated basket with carry handle, a nudge bar, rear suspension, and rear tie-down points for secure transport.

### Features

- Gas assist adjustable tiller
- Comfortable contoured high back seat with swivel and height/depth adjustments
- USB charging port
- Adjustable rear vision mirrors

- Rear adjustable suspension
- LED headlight and horn
- Width, height and angle adjustable flip-up padded armrests
- Angle adjustable backrest and headrest
- Updated basket with carry handle and rear storage pocket
- Padded loop handlebar with fingertip controls
- Variable speed settings, forward/reverse drive
- Tie-down points for secure transport

## Specifications

- **Length:** 1280mm
- **Width:** 600mm
- **Height:** 1070mm
- **Safe Climbing Angle:** 9°
- **Max Climbing Angle:** 9°
- **Maximum Load:** 136kg
- **Charger:** 8 Amp
- **Motor Size:** 0.75Hp high torque motor
- **Maximum Range:** Up to 32km
- **Weight (with battery):** 99kg

- **Turning Circle:** 1500mm