


SHOPRIDER PUMA 14HD



SKU: 888WNLLHD



PRODUCT DESCRIPTION

The Shoprider **Puma 14HD** electric powerchair has been specifically added to cater for the user who requires extra weight capacity (205kg) in a powerchair. The Puma 14HD has been ergonomically engineered to provide the robust combination of style and comfort. The Puma 14HD has the flexibility of a multi adjustable 510mm Captains Seat i.e. Seat Cushion Height, Seat Cushion Angle, Backrest Angle & Forward or Aft Adjustment on Seat Base, when combined with the powerful smooth motors, mid wheel manoeuvrability this electric wheelchair will not disappoint. The Puma's smooth operation is further enhanced by the easy steer & operate joystick controller, 360mm (14") puncture proof tyres and front and rear elastometers.

The **Puma 14HD** is suited to operation on a variety of terrains and hard surfaces, both inside and outside, the 75amp batteries fitted also ensure it can travel long distances between charges.

Features

- Large Captain's Seat
- Sleek body design
- Variable Seat height
- Controller interchangeable both arms
- Smooth operation
- Free wheel device
- Stamp supports

- Adjustable arm rests
- Option elevating leg rests
- Option of swing away leg rest

SPECIFICATIONS

- Length: 1050mm
- Width: 640mm
- Height: 1180mm
- Battery Capacity (Minimum Amp): 2 x 75 Amp
- Maximum Range: up to 25km
- Maximum Load Single User: 205kg
- Charger off board: 8 Amp
- Motor Size: 1.34 Hp
- Maximum Speed: 5 kph
- Safe Climbing Angle: 8 Degrees
- Max Climbing Angle: 9 Degrees
- Turning Radius: 500mm
- Wheel Size
 - Front Caster: 2 x 175mm
 - Middle (drive): 2 x 350mm
 - Rear Caster: 2 x 175mm