



SHOWERBUDDY SB2 TUBBUDDY

\$3,177.00



SKU: SB2

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Helping you get into a shower when the shower has a bath that may unusable. The SB2t even tilts back to assist further.

The SB2 bath transfer system from Showerbuddy is designed to allow a user to be safely transferred by a caregiver over a bathtub to shower without costly bathroom remodelling. The SB2 TubBuddy eliminates unnecessary transfers between toileting and bathing making it comfortable for both user and caregiver, while reducing costs and contributing to good user hygiene.

The SB2T bath transfer system from Showerbuddy is also designed to allow a user to be safely transferred by a caregiver over the bathtub, but with the addition of a tilt. The SB2T offers a significant reduction in costs, no need for costly bathroom remodelling. While eliminating unnecessary transfers, the SB2T with a tilt function allows for pressure relief and alleviates caregiver strain by allowing the user to slide over the bathtub without having to lift their legs. The SB2T TubBuddy Tilt can roll straight over most toilets and can also be used as a commode chair.

FEATURES

- 160kg (approx. 350lbs) weight capacity
- Rolls over most standard toilets too
- Sturdy and low maintenance structure
- Movable and removable arms
- Adjustable height
- Built in safety mechanisms

- Frame made of stainless steel and corrosion resistant aluminium.
- Safety belts and footrest
- No need for tools to assemble
- Lockable caster wheels
- **SB2t Only:** Can tilt backwards with headrest for support.

SPECIFICATIONS

- Chair Width – over caster wheels 490mm
- Chair Width – over arms 570mm
- Chair Width – inside arms 475mm – smaller users see accessory pages
- Chair Length – including footrest 845mm (745mm excluding)
- Backrest Width 460mm
- Backrest Height from Cushion to top 440mm
- Toilet Height clearance (maximum) 535mm
- Toilet Width clearance at narrowest 375mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Bridge Length (Standard) 300mm or Accessory 200mm or customized, see Accessories
- Seat Cushion Dimensions 460mm W x 430mm D
- Seat Height from Floor 575mm to 675mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Bath Tub Base unit width over feet 335mm
- Bath Tub Height from floor 535mm maximum or use RBE24 Chair Height Extender accessory
- Bath Tub Depth up to 530mm maximum or use TBE24 Bath Base Height Extender accessory
- Bath Tub width for stabilisers 440mm – 680mm
- Box Dimensions L 680mm x W 580mm x H 810mm
- Box Gross Weight 32.4kg

WILL THIS FIT YOUR BATHROOM?

[Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!](#)