

SHOWERBUDDY SB2T TUBBUDDY TILT

\$3,740.00



SKU: SB2T



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The SB2t TubBuddy Tilt bath transfer system is designed to allow a caregiver to transfer the user safely and securely over the bathtub, offering a significantly cheaper alternative to a full bathroom remodel.

The tilt function has the additional benefit of providing pressure relief as well as reducing caregiver burden by allowing the user to transfer into the bathtub without having to raise their legs.

The SB2t can also roll over most standard toilets, or function as a standard toilet commode.

FEATURES

- Rolls over most standard toilets
- Arms - removable / lockable / flip-back
- Height adjustable
- No rust aluminium frame (Grade 6061 T6)
- Removable and swing away footrests
- Safety belts - lap & chest belts plus calf straps
- Caster Wheels - 125mm x 4 locking
- Tool-free assembly

SPECIFICATIONS

- Overall Height - 108cm-118cm
- Overall Depth - 144cm (*when footrests extended and fully tilted*)
- Overall Width - 57cm

- Seat Height – 63cm-73cm
 - Seat Depth – 43cm
 - Seat Width – 47.5cm
 - Tilt Angle – 0-30 degrees (*infinite positioning*)
 - Toilet Clearance – 40cm-50
 - Weight Capacity – 250kg
-

WILL THIS FIT YOUR BATHROOM?

[Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!](#)